

### TECHNICAL SPECIFICATIONS

Stent Implant			
Stent Material	Cobalt Chromium, L-605		
Strut Thickness	81 µm		
Drug and Drug Dose	Novolimus, a Sirolimus analogue 5 µg / mm		
Polymer	Durable Methacrylate Coating, 3µm		
Crossing Profile*	.041" / 1.04 mm		
Post dilatation limit	Stent Diameter	Post-dilatation Limit per IFU	Bench Acute Testing Limit
	2.25 mm	2.75 mm	4.0 mm
	2.5 mm	3.0 mm	4.0 mm
	2.75 mm	3.25 mm	4.0 mm
	3.0 mm	3.5 mm	5.5 mm
	3.5 mm	4.0 mm	5.5 mm
	4.0 mm	4.5 mm	5.5 mm
Delivery System			
Balloon Material	Nylon Copolymer		
Working Catheter Length	140 cm		
Tip Entry Profile	0.017" , 0.43 mm		
Marker Bands	Platinum Iridium, swaged		
Catheter Shaft Diameter	Distal Outer Shaft OD 2.7F / 0.89 mm Intermediate Shaft OD 2.5F / 0.85 mm Proximal Hypotube OD 1.95F / 0.65 mm		
Nominal Pressure	9 atm		
Rated Burst Pressure	16 atm		
Guide Wire Compatibility	0.014", 0.36 mm		
Guiding Catheter Compatibility	5F (min I.D. 0.058", 1.47 mm)		

\*DESyne X2 2.5 mm diameter

# DESyne<sup>®</sup> X2

Novolimus Eluting Coronary Stent System

## COMPLIANCE CHART

PRESSURE	STENT DIAMETER					
	2.25mm	2.5mm	2.75mm	3.0mm	3.5mm	4.0mm
6 ATM		2.31	2.54	2.80	3.30	3.76
7 ATM	2.14	2.38	2.62	2.88	3.38	3.86
8 ATM	2.22	2.45	2.70	2.94	3.45	3.94
9 ATM (NOM)	2.28	2.52	2.77	3.00	3.51	4.02
10 ATM	2.35	2.58	2.84	3.06	3.56	4.09
11 ATM	2.40	2.64	2.91	3.11	3.61	4.16
12 ATM	2.45	2.70	2.97	3.15	3.66	4.22
13 ATM	2.49	2.75	3.02	3.20	3.71	4.27
14 ATM	2.53	2.80	3.06	3.24	3.75	4.32
15 ATM	2.57	2.84	3.11	3.28	3.80	4.37
16 ATM (RBP)	2.60	2.88	3.14	3.33	3.84	4.42
17 ATM	2.63	2.92	3.18	3.37	3.88	4.47
18 ATM	2.66	2.96	3.21	3.42	3.93	4.52

### FOR MORE INFORMATION OR TO ORDER THE DESyne X2 STENT SYSTEM

Contact your local sales representative or email us at [Sales@elixirmedical.com](mailto:Sales@elixirmedical.com)

Elixir Medical Corporation  
920 N. McCarthy Blvd.  
Milpitas, CA 95035 USA.  
408.636.2000  
[info@elixirmedical.com](mailto:info@elixirmedical.com)  
[www.elixirmedical.com](http://www.elixirmedical.com)

Elixir and DESyne are registered trademarks in the US and internationally for Elixir Medical Corporation.  
DESyne X2 Novolimus Eluting Coronary Stent System is CE Mark approved.  
©2024 Elixir Medical Corporation. All rights reserved.  
PMN 470 Rev F

