




HERTZ CONTACT INTRAVASCULAR LITHOTRIPSY

BREAK FREE

LithiX™ Hertz Contact Intravascular Lithotripsy Catheter

TECHNICAL SPECIFICATIONS

LithiX HC-IVL Catheter	
Available Balloon Nominal Diameters	2.0, 2.5, 3.0, 3.5, 4.0 mm
Catheter Working Length	140 cm
Catheter Design	Rapid Exchange
Guiding Catheter Compatibility	6 French (≥ 0.071 " ID or 1.8 mm ID)
Guidewire Compatibility	0.014" or 0.36 mm



Elixir Medical Corporate Headquarters
920 N. McCarthy Blvd, Milpitas, CA 95035 USA
+1.408.636.2000
info@elixirmedical.com

Elixir is a registered trademark and LithiX is a trademark in the US and internationally for Elixir Medical Corporation.

LithiX Hertz Contact Intravascular Lithotripsy Catheter is CE Mark approved.

©2026 Elixir Medical. All rights reserved.

PMN 2050 Rev C. International (OUS) Use Only.

www.ElixirMedical.com





HERTZ CONTACT INTRAVASCULAR LITHOTRIPSY

BREAK FREE

COMPLIANCE CHART

Inflation ATM	Balloon Outer Diameter (mm)				
	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm
5 [NOM]	2.06	2.57	2.99	3.59	4.03
6	2.09	2.61	3.04	3.65	4.10
7	2.12	2.65	3.09	3.70	4.16
8	2.15	2.69	3.14	3.76	4.21
9	2.18	2.72	3.18	3.80	4.26
10	2.20	2.75	3.21	3.83	4.31
11	2.22	2.77	3.24	3.86	4.35
12 [RBP]	2.24	2.79	3.27	3.88	4.40

PRODUCT ORDERING INFORMATION

LithiX Diameter				
2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm
DAR2014	DAR2514	DAR3014	DAR3514	DAR4014

FOR MORE INFORMATION OR TO ORDER THE LITHIX HC-IVL CATHETER:
Contact your local sales representative or email us at Sales@elixirmedical.com

Elixir Medical Corporate Headquarters
920 N. McCarthy Blvd, Milpitas, CA 95035 USA
+1.408.636.2000
info@elixirmedical.com

Elixir is a registered trademark and LithiX is a trademark
in the US and internationally for Elixir Medical Corporation.

LithiX Hertz Contact Intravascular Lithotripsy Catheter is CE Mark approved.
©2026 Elixir Medical. All rights reserved.
PMN 2050 Rev C. International (OUS) Use Only.

www.ElixirMedical.com

